



## Job Aid

Document G\_90\_J\_21\_A Author: T Ikanofi Authorised: W Porau 6/10/22 Review: 6/10/25

### Maldi-TOF Harsh Cleaning and Reagent Preparation Guide

#### Preparation of 80% aqueous Trifluoroacetic Acid (TFA)

80% aqueous TFA solution is used to clean MALDI-TOF target plates.

1. Transfer 400 $\mu$ L of HPLC-grade water (deionised water) into a 1.5mL eppendorf plastic tube
2. In the fume hood, carefully add 1600 $\mu$ L (1.6mL) of 99% trifluoroacetic acid (TFA) to the tube
3. Close the lid tightly on tube
4. Mix by inverting the tube 5 times
5. Label tube with TFA, date and lot number of neat TFA bottle

#### Target Cleaning Protocol

Bruker MBT target plates are reusable and will need to be cleaned regularly. To clean the targets:

1. In a petri dish, cover surface of MALDI target plate with 70% aqueous ethanol
2. Incubate at room temperature for 5 minutes
3. Remove the target and rinse it intensively under flowing tap water
4. Wipe the target (Kimwipe) intensively with 70% aqueous ethanol
5. Rinse the target plate with tap water and wipe target with a kimwipe
6. **In the fume hood**, Add 100 $\mu$ L of 80% Trifluoroacetic Acid (TFA) to the slide and wipe over all target spots on the target plate
7. Rinse the target with deionised water and dry plate with a Kimwipe
8. Let target completely dry for 15 minutes
9. Plate is now clean and ready to use

**Harsh cleaning protocol using 80% aqueous Trifluoroacetic Acid (TFA)**

