

Job Aid

Document G 90 J 10 A Author: J Ferguson Authorised: W Porau 7/04/22 Review: 7/04/24

Bacterial Isolates for Referral

For each referred isolate, complete the Inventory of stored bacterial and fungal isolates form (G_90_WS_2) and record isolate on the as a "sendaway" in the LIMS.

Fresh subcultures of isolates specified below are to be stored in glycerol broth at -20 deg C pending transport in the monthly air shipment to the central reference lab at CPHL/PMGH.

1. Any species of the Enterobacterales family (*E. coli, Klebsiella, Enterobacter* etc) isolated from any specimen site that tests as resistant to meropenem. <u>Biobank: all presumptive Carbapenemase producing Enterobacterales (CPE) – flag is meropenem resistance</u>

2. Bloodstream isolates as specified:

Use 14 day cut-off – i.e. repeat isolates of same species from the same patient not sent if < 14 days.

- Staphylococcus aureus (VRSA, MRSA or MSSA)
- Other Staphylococcus species (coagulase negative) maximum 5 isolates per month
- Enterococcal species (E. faecium or E. faecalis)
- Streptococcus pneumoniae
- Streptococcus Viridans (other alpha-haemolytic strep. species) isolated from 2 or more different blood culture collections from same patient within 2 days
- Penicillin resistant beta-haemolytic Streptococci groups A, B, C or G (PEN1 disc zone < 18mm)
- E. coli, Klebsiella, other Enterobacterales species, Salmonella (Biobank: ESBL isolates)
- Acinetobacter or Pseudomonas species (see below) Biobank: Acinetobacter baumannii complex
- Neisseria meningitidis
- Other Gram negative diplococci/coccobacilli (unidentified)

3. CSF isolates:

- Streptococcus pneumoniae
- Neisseria meningitidis
- Any other isolate obtained from a turbid CSF sample

4. Faeces

- Vibrio cholerae
- <u>Salmonella</u> or <u>Shiqella</u> species

5. Isolates from other specimens

- <u>ESBL E. coli, Klebsiella</u> species from urine if isolate is resistant to either ceftriaxone OR ceftazidime on disc testing
- Up to 10 other Gram negative isolates of any species type from urine per month
- Acinetobacter species (presumptive) isolates from blood, sterile site (i.e. pleural, peritoneal, joint fluid etc) or ET aspirate (ventilated patient – PMGH only) <u>Biobank: Acinetobacter baumannii</u> complex only
- *Pseudomonas aeruginosa* isolates from blood, sterile site or ET aspirate (ICU) <u>Biobank: meropenem-resistant isolates</u>
- Neisseria gonorrhoeae
- All suspected <u>Burkholderia pseudomallei</u> or Bacillus anthracis

Reference Lab. biobanking required for species / phenotypes that are underlined.

Reference: G 90 SOP 5 Bacterial Isolate Preparation for Storage and or Transport 11Feb22.