



National Department of Health

Title: Disc diffusion AST quality control

ID: G_90_SOP_3_A

Developed by: V Fabila, T Ikanofi, F Courtney
Reviewed by: J Ferguson, C Allen
Authorized by: W Porau
Issued: 18/11/21
Review: 2 years

Changes to the last authorized version:

Version	Date issued	Changes
A	19/11/21	N/a- New document

Certification of printed copy:

Version	
Authorised by (name)	
Signed	
Date	

NB. Printed copies of this document for local laboratory use require endorsement on the front page and recording on the Printed Controlled Document Log sheet, G_10_WS_5.

1. Purpose and Scope

ATCC QC bacterial strains are used to monitor the performance of antibiotic discs and media used for disc diffusion so that reliable antimicrobial susceptibility results are reported.

2. Principle/Clinical application

Properly stored ATCC strains exhibit predictable antimicrobial resistance properties and serve as a reliable test of the entire AST system.



Title: Disc diffusion AST quality control

ID: G_90_SOP_3_A

Revision Number: 1

Issue date: 18/11/21

Page 2 of 7

3. Responsibilities

Role	Responsibility
Quality officer or delegate	Prepares weekly QC runs and records results Highlights problems that arise without delay in the QC summary sheet Follows up on agreed corrective actions
Micro Lab manager	Reviews summary results with the QO and pathologist - discusses necessary corrective action(s) if required Signs off on results and actions before filing Tables QC records at the next Staff Meeting
Pathologist	Reviews summary QC record, signs and agrees necessary corrective actions

4. Specimen

Not applicable

5. Equipment/Materials

- stock cultures of ATCC strains renewed monthly
- agar plates
- antibiotic disc sets
- ruler for measurement

6. Procedure

Performed once a week:

Monday: subculture each ATCC organism from your current working culture as follows:

- *E. coli* ATCC 25922 –blood agar
- *S. aureus* ATCC 29213 –blood agar
- *E. faecalis* ATCC 29212 –blood agar
- *S. pneumoniae* ATCC 49619 –blood agar
- *H. influenzae* ATCC 49766 – chocolate agar



Title: Disc diffusion AST quality control

ID: G_90_SOP_3_A

Revision Number: 1

Issue date: 18/11/21

Page 3 of 7

Tuesday (at 2pm): set up disc diffusion for each combination of antibiotics using ATCC strains as follows:

Stamp	STAPH	ENC	GNR	MGRN	SSV	HAEM	STREP
Control strain	ATCC 29213	ATCC 29212	ATCC 25922	ATCC 25922	ATCC 25922	ATCC 49766	ATCC 49619
	S. aureus	E. faecalis	E. coli	E. coli	E. coli	H. influenzae	S. pneumo
Media	MH	MH	MH	MH	MH	MH-F	MH-F

Wednesday (morning):

- Read zone diameters in mm for each antibiotic/organism combination using a ruler
- Plot each zone diameter on the disc diffusion QC chart, along with the date, the disc batch number, and your initial (See Figure 2).
- If the zone diameter falls within the acceptable range (blue or green shaded areas), it is acceptable to use that disc to test patient samples.
- If the zone diameter does not fall within the acceptable range, you cannot use your disc and/or media to reliably test patient isolates. If your disc is out of range, refer to Section 8 for possible actions.
- Prepare the weekly summary QC zone sheet (Appendix for example)
- Meet with the lab manager and/or pathologist ASAP if necessary to discuss non-conforming (out of range) results



Title: Disc diffusion AST quality control

ID: G_90_SOP_3_A

Revision Number: 1

Issue date: 18/11/21

Page 5 of 7

8. Interpretation

Factors affecting zone sizes: zones too small

Observation	Probable cause	Corrective action
AMP, AUG, CRO	Disc has lost potency	Check storage conditions and package integrity.
Any zones too small	Contamination	Pick well isolate colonies when preparing inoculum. Subculture/ID inner zones to check purity and repeat if necessary
Any zones too small	Inoculum too heavy	prepare 0.5 McFarland.
Any zones too small	Too long between plate inoculation and putting discs on	Adhere to 15-15-15 minute rule
Any zones too small	Incorrect endpoint reading. E.g use of magnification	Follow EUCAST recommendations for endpoint read. i.e. in most cases measure zone diameter with unaided eye
Any zones too small	Discs not applied properly	Make sure discs are pressed firmly against the agar. Do not move a disc once applied

Zones too large

Observation	Probable cause	Corrective action
Penicillins, tetracycline, quinolones	pH of media too low	Avoid CO ₂ incubation for MHA (clear plates) – lowers pH of media
Any zones too large	Inoculum too light	prepare inoculum to 0.5 McFarland standard. Make sure your lawn of growth is confluent with no spaces
Any zones too large	Incorrect endpoint read	Take into account double zones, colonies within the zone. Read at point of complete inhibition with fuzzy zone edges
Any zones too large	Inoculum source plate too old and contains too many non-viable cells	Use overnight growth – subculture if necessary. When reading sensitivity plates, should have a confluent lawn of growth with no spaces.

9. Safety

As per standard precautions.

10. Quality Control

Not applicable.



Title: Disc diffusion AST quality control

ID: G_90_SOP_3_A

Revision Number: 1

Issue date: 18/11/21

Page 6 of 7

11. Related documents

For access, refer to <https://path-png.org/microbiology-sops-fleming-fund/>

Weekly AST QC summary	G_90_WS_3
Weekly AST QC Worksheets	G_90_WS_4
Antibiotic disc susceptibility testing	G_90_SOP_6
Maintenance of cultures used for quality control testing	G_90_SOP_2

12. Reference

European Committee on Antimicrobial Susceptibility Testing (EUCAST) QC tables: https://www.eucast.org/ast_of_bacteria/quality_control/



Title: Disc diffusion AST quality control

ID: G_90_SOP_3_A

Revision Number: 1

Issue date: 18/11/21

Page 7 of 7

Appendix: Weekly Disc AST QC summary results (example)

		Worksheet (form)	
Document G_90_WS_3_A Author: F Courtney Authorised: W Porau 04/08/2021 Review: 04/08/21			
Weekly AST QC Summary			

Week commencing: 15/7/20

Prepared by: V Fabila Signed: [Signature]

Summary assessment: _____

(Tick cells where zone size is within control limits; otherwise indicate zone size readings with expected results in brackets)

Stamp	STAPH	ENC	GNR	MRGN	SSV	HAEM	STREP
Control strain	ATCC 29213	ATCC 29212	ATCC 25922	ATCC 25922	ATCC 25922	ATCC 49766	ATCC 49619
	S. aureus	E. faecalis	E. coli	E. coli	E. coli	H. influenzae	S. pneumo
Media	MH	MH	MH	MH	MH	MH-F	MH-F
P1	✓					✓	✓
OX1							✓
E15	✓						✓
TE30	✓					✓	✓
SXT25	✓		✓		✓	✓	✓
C30	✓			✓	✓	✓	✓
CRO30			✓		✓	✓	
FOX30	✓						
AMP2		✓					
F100		✓	✓				
VA5		✓					
AMC30			✓				
CN10			17 (19-26)				
CIP5			✓				
TOB10				✓			
MEM10				✓			
CAZ10				✓			
TZP30/6				✓			
AK30				✓			
AMP10					✓		
PEF5					✓		
AZ15					✓		

Action recommended (lab manager):

Repeat gent. with new batch.

Action required (supervisor):

Supervisor name: Dr Joseph Signed: [Signature] Date: 20/7/20